

1 Choose sensor site by age

(If worn on sites not listed below, sensor may not work as expected. Talk to your healthcare provider about the best site for you.)

- a. Age 7 years and older: Use back of upper arm or abdomen.
- b. Age 2 years (24 months) – 6 years: Use back of upper arm, upper buttocks, or abdomen.

Avoid areas:

- With loose skin or without enough fat to avoid muscles and bones
- That get bumped, pushed, or you lie on while sleeping
- Within 8 centimeters of infusion or injection site
- Near waistband or with irritations, scarring, tattoos, or lots of hair

To keep the sensor attached see dexcom.com/faqs

2 Clean and dry site

- a. Wash hands with soap and water. Dry.
- b. Rub site with an alcohol wipe to get rid of all oils.
- c. Let dry completely so sensor will stick.

3 Unscrew cap

- Don't touch inside applicator.
- Don't use if damaged or previously opened.

4 Use applicator to insert sensor

Relax any muscles near site.

- ① Press and hold applicator firmly against skin until clear safety guard is pushed in.
- ② Press the button while the guard is pushed in.

Tip: Pressing applicator firmly against the skin unlocks the button.

5 Remove applicator

- a. Save applicator because you'll need its pairing code.
- b. Screw cap back on applicator.

6 Rub patch

- a. Rub firmly around patch 3 times.
- b. Gently press on top of sensor for 10 seconds.

7 Apply overpatch (when needed)

Apply overpatch when needed to help keep the sensor on your body. Find the overpatch bundled with the instructions in the sensor box.

- a. Carefully pull off both clear liners, one at a time. Don't touch white adhesive area.
- b. Use colored tab to place overpatch around sensor.
- c. Rub around overpatch.
- d. Use tab to peel off colored liner.
- e. Rub around overpatch again.

8 Continue setup

- Go back to the app or receiver to finish setup.
- For tips on keeping the sensor attached, go to dexcom.com/faqs.

Adhesive patch care

- The longer you keep the patch dry and sweat-free in the first 12 hours, the longer it may stick to your skin.
- When it gets wet, gently pat it dry as soon as you can.
- If it peels off your skin, trim the peeled parts and put on medical tape.

Contraindications

No MRI/CT/diathermy — MR unsafe: Don't wear

Interpretation of Dexcom ONE+ results should be based on glucose trends and several sequential sensor readings over time. Dexcom ONE+ also in the detection of episodes of hyperglycemia and hypoglycemia, facilitating both acute and long-term therapy adjustments.

Dexcom ONE+ is intended for use by patients at and in healthcare facilities.

Warnings

Before using this product, read the instructions the Dexcom ONE+ User Guide, available at dexcom.com/downloadsandguides. Don't ignore low or high symptoms: Use your meter to make treatment decisions when your readings don't match your low or high symptom needed, seek immediate medical attention.

No number or no arrow, no CGM treatment d

Use your Bg meter to make treatment decisions your Dexcom ONE+ doesn't show both a number and trend arrow. Also use your Bg meter during 30-minute sensor warmup period.

Don't use if you are on dialysis or critically ill. Dexcom ONE+ performance hasn't been evaluated in these populations and sensor readings may be inaccurate.

Sensor wire breaks off: Don't ignore broken or detached sensor wires. If this happens, please contact Dexcom Technical Support at dexcom.com or your local Dexcom representative.

Packaging Symbols

Symbols on the Dexcom ONE+ packaging show may not have meaning in your region and are dexcom.com for symbols not defined here.

Symbol Definition

IPXX

UK responsible person

UK CA

UK marking of conformity

UK REP

R-NZ

New Zealand radio compliance